



Surgicoll® collagen haemostat

An absorbable, highly purified collagen haemostat sponge for use in surgical procedures

MBP Surgicoll® provides quick and reliable haemostatic action, activating the clotting cascade. It is a white, sponge-like, easily mouldable and insoluble fleece. It conforms to the wound surface and has excellent adhesion properties. Surgicoll® contains no formaldehyde or other chemicals and is simply collagen and a small percentage of water. It can be applied to the wound surface with dry instruments and haemostasis is achieved in one to four minutes. It may be cut to shape and size with sterile scissors. Absorbency is between two and four weeks depending on the amount used. Surgicoll® has excellent tissue adherence and may be sutured in place. Shelf boxes are designed to conveniently fit suture racks for convenient access at any time during surgery

Product Ordering Information:

Cat. No.	Size (cm)	Quantity Per box
AV0050	1 x 1	30
AV0051	2 x 2	12
AV0052	5 x 5	10
AV0053	5 x 10	10
AV0054	10 x 10	5

Surgicoll is manufactured by Medical Biomaterial Products GmbH

PROBRO04R1