



Nuros launches AndraSnare

Advanced design sets new performance standards in retrieval:

- Highly visible, gold plated, atraumatic loop, resists collapse
- Platinum inner marker for accurate positioning
- Reinforced catheter to hub joint avoids kinking
- Minimal back-bleeding from convenient insertion tool
- Large, positive torque device and almost 1:1 torque
- Wide range, including 'micro' sizes, from 2 to 35mm



AndraSnare



LoopMaster

LoopMaster-Sochman Snare


LoopMaster cleverly overcomes the challenge of snaring when an open end is difficult to access:

- Locate the snare loop over the object laterally
- Advance the integral wire around the opposite side of the object and into the snare loop
- Pull back the snare loop trapping the object firmly


M.A.R.S. Multilayer Flow Modulators

The aneurysm repair system that achieves physiological exclusion whilst maintaining flow to collaterals and side branches:

- Diameters 28 to 40mm now available for TAAA, CE marked.
- 54 patient, multicentre peripheral/visceral study: 12 month follow up published, reporting:
- No sac reperfusion, no aneurysm enlargement
- No collateral vessel obstruction,



M.A.R.S.



X-Cath and HQS Sheath

X-Cath Multi-Port Adapter

Conveniently provides multiple entry via single sheath access for complex endovascular procedures:

- Accommodates up to 3 devices each 9F max.
- Connects to the HQS introducer sheath
- Ensures minimal blood loss at all stages
- HQS sheath available in 18 to 26F ID.
- Charing Cross Innovation prize winners